





qube®

Qube[®] facilitates at-a-glance monitoring with an ultra-bright display featuring large numerics. Multicolor, severity-based lights illuminate on both the front and back of the monitor for enhanced alarm surveillance.



qube®



The Qube compact monitor delivers both simplicity and style.

INTRODUCING QUBE

Small and lightweight with a battery life that goes the distance, Qube offers portability, accessibility and heightened connectivity—perfect for emergency-, general-, intensive- and post-anesthesia care units. Caregivers can access current, critical patient information— from bedside throughout transport—to provide an informed level of patient care.

VERSATILE

Qube is mountable anywhere – wall, roll-stand, bed rail, anesthesia system – or can quickly be detached via a quick-release feature for immediate mobility. Long battery life and dual battery capability keep Qube on the go everywhere your patient goes. Impact resistant, chemical resistant, and protected from liquid ingress, Qube endures even with heavy use. Full functionality in a compact package.





SPACELABS
HEALTHCARE

INSIGHTFUL

Data from a variety of sources informs and supports clinical determinations at the bedside or remotely. Graphical assessment tools and the ability to access up to 96 hours of physiologic data empowers clinicians and helps facilitate decision making.

CONFIDENT

Clinicians can customize the screen view to quickly access a multitude of graphical and tabular trend presentations based on individual patient needs and procedural guidelines. Interactive view, control, and recording for parameters, trends, and calculations can be engaged instantaneously, linking clinicians with a range of real-time patient physiologic data. Qube enables point-of-care decisions that provide ongoing reassurance and confidence to patients and their loved ones.


ENLIGHTENED

Qube facilitates at-a-glance monitoring with an ultra-bright display featuring large numerics. Multicolor, severity-based lights illuminate on both the front and back of the monitor for enhanced alarm surveillance.

ATTRACTIVE

The sleek, contemporary details of Qube, including ergonomic user controls and a comfortable grip handle, create a positive impact for caregiver and patient. Art meets technology.

The small, compact form of the capnography pod makes it ideal for use with Qube. Designed to attach seamlessly onto the rear of the monitor, it's ready in a snap.



CONNECTED AND COMPATIBLE

With industry-standard communication protocols Qube provides inter-monitor communication, enhanced alarm surveillance, and broader productivity through both wired and wireless channels. The Qube is also designed for compatibility with the Spacelabs Command Modules, Capnography Pod, and Exergen Temporal Artery Thermometer to provide a versatile, full range of measurement choices. With Spacelabs parameter modules, you can easily select from a variety of configurations to suit the monitoring needs of specific patients or care units in your hospital.

CAPNO POD

The new Capno Pod measures end-tidal carbon dioxide (EtCO₂), minimum CO₂, and respiration rate to help determine your patient's ventilatory, circulatory, and metabolic status. Designed especially for the Qube monitor, it is suitable for use on patients of all ages.



SHAPING THE FUTURE OF PATIENT MONITORING

A Quick Look at the Qube Compact Patient Monitor. Portability, Easy to Use, Heightened Connectivity.

- 3, 5, or Diagnostic 12 lead ECG, Resp, SPO2, NIBP, Temp
- 4 invasive pressures and Cardiac Output
- Up to 6 waveforms, 18 parameters
- Optional Secondary 15" (38 cm) or optional 19" (48.3 cm) interactive touchscreen display
- Standard 8 hour battery using 2 hot swappable batteries
- 24 hour Patient Data Shuttle
- One button, quick release Docking Station
- Wireless connectivity
- Internal strip recorder
- Optional Exergen USP Spot Check Thermometer
- Optional Capnography
- LifeSync compatible

Features may vary depending on purchased configuration.

